

## Dioxopromethazine hydrochloride

Catalog No: #S3723


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Dioxopromethazine hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Histamine Receptor inhibitor
Calculated MW	348.847
Formulation	C17H20N2O2SHCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: CC(CN(C)C)CC1=CC=CC=C1S(=O)(=O)C2=CC=CC=C2  
Pathways: GPCR/G Protein; Immunology/Inflammation; Neuroscience  
Receptor: H1 receptor  
Boiling pt: Melting pt: Solubility: Soluble in DMSO  
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. He LL, et al. J Fluoresc. 2011 Sep;21(5):1847-56.

Note: This product is for in vitro research use only