

Rolapitant hydrochloride

Catalog No: #S3724


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rolapitant hydrochloride
Brief Description	Inhibitors
Purification	99.79%
Target Name	Neurokinin 1 antagonist
Calculated MW	536.94
Formulation	C25H27ClF6N2O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous

system SMILES: O=C1N[C@@]2(CN[C@@](C3=CC=CC=C3)(CO[C@@H](C4=CC(C(F)F)=CC(C(F)F)=C4)C)CC2)CC1.[H]Cl
 Pathways: Neuroscience
 Receptor: NK1
 Boiling pt: Melting pt: Solubility: Soluble in DMSO
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Rapoport B, Schwartzberg L, Chasen M, Powers D, Arora S, Navari R, Schnadig I. Eur J Cancer. 2016 Apr;57:23-30.

Note: This product is for in vitro research use only