

K-Ras(G12C) Inhibitor 6

Catalog No: #S3725


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	K-Ras(G12C) Inhibitor 6
Brief Description	Inhibitors
Purification	98.00%
Target Name	Raf inhibitor
Calculated MW	405.334
Formulation	C17H22Cl2N2O3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:CCCCNC(=O)[C@H](C)C[C@@H]([C@H](CC(C)C)NC(=O)[C@H](C)NC(=O)[C@H](CC(C)C)NC(=O)C)OPathways:MAPKRecept
 or:K-RasBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and
 shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

- Hunter JC, Gurbani D, Ficarro SB, Carrasco MA, Lim SM, Choi HG, Xie T, Marto JA, Chen Z, Gray NS, Westover KD. Proc Natl Acad Sci U S A. 2014 Jun 17;111(24):8895-900.

Note: This product is for in vitro research use only