

Methyl syringate


Catalog No: #S3727

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Methyl syringate
Brief Description	Inhibitors
Purification	99.00%
Target Name	TRP/TRPV Channel inhibitor
Calculated MW	212.2
Formulation	C10H12O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC(=O)C(=C(C=C1OC)OC)C(=O)OC
Pathways: Membrane transporter/Ion channel
Receptor: hTRPA1
Boiling pt: Melting pt: Solubility: Soluble in Chloroform, Dichloromethane, Ethyl Acetate, Acetone, DMSO
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Kim MJ, et al. PLoS One. 2013 Aug 21;8(8):e71603.

Note: This product is for in vitro research use only