

Methyl gallate


Catalog No: #S3730

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Methyl gallate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Reverse Transcriptase inhibitor
Calculated MW	184.14
Formulation	C ₈ H ₈ O ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: COC(=O)C1=CC(=C(C(=C1)O)O)O
Pathways: Microbiology & virology
Receptor: Reverse Transcriptase
Boiling pt: Melting pt: Solubility: Soluble in DMSO
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Wang CR, et al. Environ Toxicol Pharmacol. 2014 Mar;37(2):626-37.

Note: This product is for in vitro research use only