

3-Phenylpropionic acid


Catalog No: #S3731

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	3-Phenylpropionic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	150.17
Formulation	C9H10O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C1=CC=C(C=C1)CCC(=O)O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Soluble in DMSO
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Johannes G, et al. Tetrahedron Letters, 1998, 39(45): 8329-8332

Note: This product is for in vitro research use only