

Sebacic acid


Catalog No: #S3732

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sebacic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	HIV Protease inhibitor
Calculated MW	202.25
Formulation	C10H18O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:C(CCCCC(=O)O)CCCC(=O)OPathways:Microbiology&virology; Proteases/ProteasomeReceptor:HIVBoling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Clemente Bretti, et al. Eng. Data. 51 (5): 1660C1667.

Note: This product is for in vitro research use only