

## Methyl trans-p-Coumarat

Catalog No: #S3733


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Methyl trans-p-Coumarat
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	178.18
Formulation	C10H10O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:COC(=O)/C=C/C1=CC=C(C=C1)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Chloroform, DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1.Tan EM, et al. J Phys Chem B. 2013 May 2;117(17):4798-805.

Note: This product is for in vitro research use only