

## 3',4'-Dihydroxyacetophenone

Catalog No: #S3734


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	3',4'-Dihydroxyacetophenone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Platelet aggregation inhibitor
Calculated MW	152.14
Formulation	C <sub>8</sub> H <sub>8</sub> O <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: CC(=O)C1=CC(=C(C=C1)O)O  
Pathways: Others  
Receptor: Platelet aggregation  
Boiling pt: Melting pt: Solubility: Soluble in Chloroform, DMSO  
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Sun YP, et al. Zhonghua Fu Chan Ke Za Zhi. 1992 Jan;27(1):9-11, 57.

Note: This product is for in vitro research use only