

Tetrahydro Curcumin

Catalog No: #S3735


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tetrahydro Curcumin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	372.41
Formulation	C21H24O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: COC1=C(C=CC(=C1)CCC(=O)CC(=O)CCC2=CC(=C(C=C2)O)OC)O
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 3.6 mg/ML
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Kang N, et al. Food Chem Toxicol. 2014 May;67:193-200.

Note: This product is for in vitro research use only