

## Olivetol


Catalog No: #S3736

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Olivetol
Brief Description	Inhibitors
Purification	99.00%
Target Name	IκB/IKK inhibitor
Calculated MW	180.24
Formulation	C11H16O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: CCCCCC1=CC(=CC(=C1)O)O  
Pathways: NF-κB Receptor; IκB  
Boiling pt: Melting pt: Solubility: Soluble in Dichloromethane, Ethyl Acetate  
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Sarah K. Oettl, et al. PLoS One. 2013; 8(10): e76929.

Note: This product is for in vitro research use only