

3-(4-Hydroxyphenyl)-1-propanol

Catalog No: #S3738


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	3-(4-Hydroxyphenyl)-1-propanol
Brief Description	Inhibitors
Purification	98.00%
Target Name	LTR inhibitor
Calculated MW	152.19
Formulation	C ₉ H ₁₂ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
SMILES: C1=CC(=CC=C1CCCO)O
Pathways: Immunology/Inflammation
Receptor: LTR
Boiling pt: Melting pt: Solubility: Soluble in Chloroform, Dichloromethane, Ethyl Acetate, Acetone, DMSO
Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Crowley JJ, et al. Neurotoxicology. 2013 Dec;39:153-7.

Note: This product is for in vitro research use only