

## 2,6-Dimethoxy-1,4-benzoquinone


Catalog No: #S3739

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	2,6-Dimethoxy-1,4-benzoquinone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	168.148
Formulation	C <sub>8</sub> H <sub>8</sub> O <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: O=C(C(OC)=C1)C(OC)=CC1=O

Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Methanol, acetone Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

N/A

Note: This product is for in vitro research use only