

CtIP antibody

Catalog No: #22632

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CtIP antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 452 and 747 of CtIP
Target Name	CtIP
Accession No.	Swiss-Prot:Q99708Gene ID:5932
Uniprot	Q99708
GeneID	5932;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 102kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:500

Images



Sample (30 ug of whole cell lysate)
 A: Hep G2
 7.5% SDS PAGE
 Primary antibody diluted at 1: 5000

Background

The protein encoded by this gene is a ubiquitously expressed nuclear protein. It is found among several proteins that bind directly to retinoblastoma

protein, which regulates cell proliferation. This protein complexes with transcriptional co-repressor CTBP. It is also associated with BRCA1 and is thought to modulate the functions of BRCA1 in transcriptional regulation, DNA repair, and/or cell cycle checkpoint control. It is suggested that this gene may itself be a tumor suppressor acting in the same pathway as BRCA1. Three transcript variants encoding two different isoforms have been found for this gene. More transcript variants exist, but their full-length natures have not been determined. [provided by RefSeq]

Note: This product is for in vitro research use only