

## 4-hydroxyphenylacetic acid

Catalog No: #S3743


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	4-hydroxyphenylacetic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	152.15
Formulation	C8H8O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:OC(=O)Cc1ccc(cc1)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSOEthanolWaterAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only