

Veratric acid


Catalog No: #S3744

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Veratric acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	182.17
Formulation	C9H10O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:COC1=C(C=C(C=C1)C(=O)O)OC

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:ChoroformDMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Lee KE, et al. A study of facial wrinkles improvement effect of veratric acid from cauliflower mushroom through photo-protective mechanisms against UVB irradiation. Arch Dermatol Res. 2016 Apr;308(3):183-92.

Note: This product is for in vitro research use only