

## 3,4-Dimethoxyphenol


Catalog No: #S3745

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	3,4-Dimethoxyphenol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	154.16
Formulation	C8H10O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: C1=CC(=CC=C1CC(=O)O)O Pathways: Others Receptor: Others Boiling pt: Melting

pt: Solubility: Chloroform DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

N/A

Note: This product is for in vitro research use only