

# triacontan-1-ol

Catalog No: #S3764


Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	triacontan-1-ol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	438.825
Formulation	C30H62O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OthersSMILES:CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCOPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in

Chloroform, Dichloromethane, Ethyl Acetate, DMSO, Acetone, etc.Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only