

Armepavine

Catalog No: #S3775


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Armepavine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others; Phospholipase
Calculated MW	313.39
Formulation	C19H23NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: CN1CCC2=CC(=C(C=C2[C@H]1CC3=CC=C(C=C3)O)OC)OC Pathways: Metabolism; Others Receptor: Itk; PLCy Boling pt; Melting pt; Solubility: 10 mM in DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ka SM, et al. (S)-armepavine from Chinese medicine improves experimental autoimmune crescentic glomerulonephritis. *Rheumatology (Oxford)*. 2010 Oct;49(10):1840-51.

Note: This product is for in vitro research use only