

## Pueraria glycoside

Catalog No: #S3778


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pueraria glycoside
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	432.38
Formulation	C21H20O10
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C1=CC(=C(C=C1C2=COC3=C(C2=O)C=CC(=C3[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Jin SE,et al. Anti-inflammatory and antioxidant activities of constituents isolated from Pueraria lobata roots.Arch Pharm Res. 2012 May;35(5):823-37.

Note: This product is for in vitro research use only