

Oroxin B

Catalog No: #S3780


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Oroxin B
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	594.52
Formulation	C27H30O15
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C1=CC=C(C=C1)C2=CC(=O)C3=C(C(=C(C=C3O2)O[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO[C@H]5[C@@H]([C@H]([C@@H]([C@H](O5)CO)O)O)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

- Li DQ,et al. Discovery of xanthine oxidase inhibitors from a complex mixture using an online, restricted-access material coupled with column-switching liquid chromatography with a diode-array detection system. Anal Bioanal Chem. 2014 Mar;406(7):1975-84.

Note: This product is for in vitro research use only