

Fraxin

Catalog No: #S3783


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fraxin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	370.31
Formulation	C16H18O10
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:COC1=C(C(=C2C(=C1)C=CC(=O)O2)OC3C(C(C(C(O3)CO)O)O)O)OPathways:OthersReceptor:mURAT1 or mGLUT9Boiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Whang WK,ET AL. Natural compounds,fraxin and chemicals structurally related to fraxin protect cells from oxidative stress. Exp Mol Med. 2005 Oct 31;37(5):436-46.

Note: This product is for in vitro research use only