

Griffonin

Catalog No: #S3784


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Griffonin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	329.3
Formulation	C14H19NO8
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C1=C/C(=C/C#N)/[C@@H]([C@H]([C@@H]1O)O)O[C@H]2[C@@H]([C@H]([C@@H]([C@H](O2)CO)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Tih AE,et al. Minor biflavonoids from *Lophira alata* leaves. J Nat Prod. 2006 Aug;69(8):1206-8.

Note: This product is for in vitro research use only