

## Ajugol

Catalog No: #S3792


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ajugol
Brief Description	Inhibitors
Purification	99.15%
Target Name	Others
Calculated MW	348.4
Formulation	C15H24O9
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C[C@@]1(C[C@H]([C@H]2[C@@H]1[C@@H](OC=C2)O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 3.7 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1.Elusiyan CA, et al. Distribution of iridoid glucosides and anti-oxidant compounds in *Spathodea campanulata* parts. *Afr J Tradit Complement Altern Med.* 2011;8(1):27-33.

Note: This product is for in vitro research use only