

## Corilagin

Catalog No: #S3795


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Corilagin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	634.45
Formulation	C27H22O18
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Others SMILES: C1[C@@H]2[C@H]([C@@H]([C@H]([C@@H](O2)OC(=O)C3=CC(=C(C(=C3)O)O)O)OC(=O)C4=CC(=C(C(=C4C5=C(C(=C(C=C5C(=O)O1)O)O)O)O)O)O

Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in

Pyridine, Methanol, Ethanol Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Wang Z et al. Corilagin attenuates aerosol bleomycin-induced experimental lung injury.

Note: This product is for in vitro research use only