

## ISORHAMNETIN 3-O-NEOHESPEROSIDE

Catalog No: #S3799


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	ISORHAMNETIN 3-O-NEOHESPEROSIDE
Brief Description	Inhibitors
Purification	99.64%
Target Name	Others
Calculated MW	624.54
Formulation	C28H32O16
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)O[C@@H]2[C@H]([C@@H]([C@H](O[C@H]2OC3=C(OC4=CC(=CC(=C4C3=O)O)O)C5=CC(=C(C=C5)O)OC)CO)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

- Li CY,et al. Measurement of hydroxysafflor yellow A in human urine by liquid chromatography-tandem mass spectrometry.J Chromatogr B Analyt Technol Biomed Life Sci.

Note: This product is for in vitro research use only