

Byakangelicin

Catalog No: #S3801


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Byakangelicin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	334.32
Formulation	C17H18O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:CC(C)([C@@H](COC1=C2C(=C(C3=C1OC(=O)C=C3)OC)C=CO2)O)Pathways:OthersReceptor:OthersBoiling pt:Melting

pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Yang J, et al. Byakangelicin induces cytochrome P450 3A4 expression via transactivation of pregnane X receptors in human hepatocytes. Br J Pharmacol. 2011 Jan;162(2):441-51.

Note: This product is for in vitro research use only