

Neomangiferin

Catalog No: #S3804


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neomangiferin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	584.483
Formulation	C25H28O16
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C1=C2C(=CC(=C1O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)OC4=CC(=C(C(=C4C2=O)O)[C@H]5[C@@H]([C@H]([C@@H]([C@H](O5)CO)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Zhou C,et al. Beneficial effects of neomangiferin on high fat diet-induced nonalcoholic fatty liver disease in rats.

Int Immunopharmacol. 2015 Mar;25(1):218-28.

Note: This product is for in vitro research use only