

Icariside I

Catalog No: #S3808


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Icariside I
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	530.526
Formulation	C27H30O11
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:CC(=CCC1=C(C=C(C2=C1OC(=C(C2=O)O)C3=CC=C(C=C3)OC)O)O[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)O)

CPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ming LG,et al. Functions and action mechanisms of flavonoids genistein and icariin in regulating bone remodeling. J Cell Physiol. 2013 Mar;228(3):513-21

Note: This product is for in vitro research use only