

Darutigenol

Catalog No: #S3817


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Darutigenol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	322.489
Formulation	C20H34O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular

systemSMILES:C[C@@]1(CC[C@@H]2C(=C1)CC[C@H]3[C@@]2(CC[C@H](C3(C)C)O)C)[C@H](CO)O

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Chen FF, et al. Molecularly imprinted polymer for the specific solid-phase extraction of kireinol from Siegesbeckia pubescens herbal extract. Talanta. 2012 Jan 30;89:505-12.

Note: This product is for in vitro research use only