

Acetylursolic acid

Catalog No: #S3819


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Acetylursolic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR inhibitor
Calculated MW	498.748
Formulation	C32H50O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: CC(C)C1CC2(C)C(C)C(=O)C3C(=O)C(C)C(C)C4(C)C(=O)C(C)C(C)C5(C)C(=O)C(C)C(C)C6(C)C(=O)C(C)C(C)C7(C)C(=O)C(C)C(C)C8(C)C(=O)C(C)C(C)C9(C)C(=O)C(C)C(C)C10(C)C(=O)C(C)C(C)C11(C)C(=O)C(C)C(C)C12(C)C(=O)C(C)C(C)C13(C)C(=O)C(C)C(C)C14(C)C(=O)C(C)C(C)C15(C)C(=O)C(C)C(C)C16(C)C(=O)C(C)C(C)C17(C)C(=O)C(C)C(C)C18(C)C(=O)C(C)C(C)C19(C)C(=O)C(C)C(C)C20(C)C(=O)C(C)C(C)C21(C)C(=O)C(C)C(C)C22(C)C(=O)C(C)C(C)C23(C)C(=O)C(C)C(C)C24(C)C(=O)C(C)C(C)C25(C)C(=O)C(C)C(C)C26(C)C(=O)C(C)C(C)C27(C)C(=O)C(C)C(C)C28(C)C(=O)C(C)C(C)C29(C)C(=O)C(C)C(C)C30(C)C(=O)C(C)C(C)C31(C)C(=O)C(C)C(C)C32(C)C(=O)C(C)C(C)C

Pathways: Neuroscience Receptor: AChE

Boling pt: Melting pt: Solubility: Soluble in Chloroform, Dichloromethane, Ethyl Acetate, DMSO

Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only