

Iridin

Catalog No: #S3820


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Iridin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	522.459
Formulation	C24H26O13
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:COC1=CC(=CC(=C1OC)O)C2=COC3=CC(=C(C(=C3C2=O)O)OC)O[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)OPa
 thways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the
 tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only