

Dihydroresveratrol

Catalog No: #S3840


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dihydroresveratrol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	230.3
Formulation	?C14H14O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:C1=CC(=CC=C1CCC2=CC(=CC(=C2)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Jung, C.M.,Heinze, T.M.,Schnackenberg, L.K.,?et al.?Interaction of dietary resveratrol with animal-associated bacteria.?FEMS Microbiology Letters. 2009, 297(2),?266-273.

Note: This product is for in vitro research use only