

Bavachin

Catalog No: #S3846


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bavachin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Estrogen/progestogen Receptor activator
Calculated MW	324.376
Formulation	C20H20O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:CC(=CCC1=C(C=C2C(=C1)C(=O)C[C@H](O2)C3=CC=C(C=C3)O)O)CPathways:Endocrinology/HormonesReceptor:ERα; ERβBoiling pt:Melting pt:Solubility:Soluble in Chloroform, Dichloromethane, Ethyl Acetate, DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Park J1,et al. Activation of Estrogen Receptor by Bavachin from Psoralea corylifolia. Biomol Ther (Seoul). 2012 Mar;20(2):183-8.

Note: This product is for in vitro research use only