

Erianin

Catalog No: #S3864


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Erianin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Bcr-Abl inhibitor
Calculated MW	318.369
Formulation	C18H22O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:COC1=C(C=C(C=C1)CCC2=CC(=C(C(=C2)OC)OC)OC)OPathways:Angiogenesis; Cytoskeletal SignalingReceptor:Bcl-2Boling pt:Melting pt:Solubility:Soluble in Chloroform, Dichloromethane, Ethyl Acetate, DMSOAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lam F, et al. ZJU-6, a novel derivative of Erianin, shows potent anti-tubulin polymerisation and anti-angiogenic activities.

Invest New Drugs. 2012 Oct;30(5):1899-907.

Note: This product is for in vitro research use only