

Anwuligan

Catalog No: #S3866


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anwuligan
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	328.408
Formulation	C ₂₀ H ₂₄ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others SMILES: C[C@H](CC1=CC2=C(C=C1)OCO2)[C@@H](C)CC3=CC(=C(C=C3)O)OC Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: Soluble in Chloroform, Dichloromethane, Ethyl Acetate, DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Yang X, et al. [Selection of mobile phases for the determination of gamma-schisandrin and multi-active constituents in Schisandra chinensis and its preparations by high performance liquid chromatography].

Note: This product is for in vitro research use only