

Alpinetin

Catalog No: #S3867


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alpinetin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Bcr-Abl inhibitor
Calculated MW	270.284
Formulation	C16H14O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:COC1=C2C(=O)C[C@H](OC2=CC(=C1)O)C3=CC=CC=C3Pathways:Angiogenesis; Cytoskeletal SignalingReceptor:Bcl-2Boling pt:Melting pt:Solubility:Soluble in Chloroform, Dichloromethane, Ethyl Acetate, DMSO, AcetoneAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Guan S,et al. Immunosuppressive activity of alpinetin on activation and cytokines secretion of murine T lymphocytes. Immunopharmacol Immunotoxicol. 2014 Aug;36(4):290-6.

Note: This product is for in vitro research use only