

## Skimmin

Catalog No: #S3869


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Skimmin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	324.285
Formulation	C15H16O8
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C1=CC(=CC2=C1C=CC(=O)O2)O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)OPathways:OthersReceptor:OthersBinding pt:Melting pt:Solubility:Pyridine, Methanol, EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Zhang S, et al. Skimmin, a coumarin, suppresses the streptozotocin-induced diabetic nephropathy in wistar rats. Eur J Pharmacol. 2012 Oct 5;692(1-3):78-83.

Note: This product is for in vitro research use only