

## Nodakenin

Catalog No: #S3872


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Nodakenin
Brief Description	Inhibitors
Purification	99.00%
Target Name	AChR inhibitor
Calculated MW	408.403
Formulation	C <sub>20</sub> H <sub>24</sub> O <sub>9</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Others SMILES: CC(C)([C@H]1CC2=C(O1)C=C3C(=C2)C=CC(=O)O3)O[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)O Pathways: Neuroscience Receptor: AChE Boiling pt: Melting pt: Solubility: Soluble in Pyridine, Methanol, Ethanol Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Park SJ, et al. Effect of nodakenin on atopic dermatitis-like skin lesions. *Biosci Biotechnol Biochem.* 2014;78(9):1568-71.

Note: This product is for in vitro research use only