

## Diosgenin glucoside

Catalog No: #S3874


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Diosgenin glucoside
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	576.771
Formulation	C33H52O8
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C[C@@H]1CC[C@@]2([C@H]([C@H]3[C@@H](O2)C[C@@H]4[C@@]3(CC[C@H]5[C@H]4CC=C6[C@@]5(CC[C@@H](C6)O[C@H]7[C@@H]([C@H]([C@@H]([C@H](O7)CO)O)O)C)C)OC1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1.Malinow MR,et al. Effects of synthetic glycosides on steroid balance in *Macaca fascicularis*. *J Lipid Res.* 1987 Jan;28(1):1-9.

Note: This product is for in vitro research use only