

Loureirin B

Catalog No: #S3876


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Loureirin B
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	316.353
Formulation	C18H20O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism/Nervous

systemSMILES:COC1=CC(=C(C(=C1)OC)CCC(=O)C2=CC=C(C=C2)O)OCP
 Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in DMSO
 Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Bai X, et al. Loureirin B inhibits fibroblast proliferation and extracellular matrix deposition in hypertrophic scar via TGF-β/Smad pathway. Exp Dermatol. 2015 May;24(5):355-60.

Note: This product is for in vitro research use only