

VERBENALIN

Catalog No: #S3890


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	VERBENALIN
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	388.369
Formulation	C17H24O10
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C[C@H]1CC(=O)[C@H]2[C@@H]1[C@@H](OC=C2C(=O)OC)[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Kang Z1,et al. Cornin induces angiogenesis through PI3K-Akt-eNOS-VEGF signaling pathway. Food Chem Toxicol. 2013 Aug;58:340-6.

Note: This product is for in vitro research use only