

## Chonglou Saponin II

Catalog No: #S3894


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Chonglou Saponin II
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	855.028
Formulation	C44H70O16
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:CC1CCC2(C(C3C(O2)CC4C3(CCC5C4CC=C6C5(CCC(C6)OC7C(C(C(C(O7)CO)OC8C(C(C(O8)CO)O)O)OC9C(C(C(C(O9)C)O)O)O)C)C)OC1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Yin X,et al. Simultaneous determination and pharmacokinetic study of polyphyllin I, polyphyllin II, polyphyllin VI and polyphyllin VII in beagle dog plasma after oral administration of Rhizoma Paridis extracts by LC-MS-MS. Biomed Chromatogr. 2013 Mar;27(3):343-8.

Note: This product is for in vitro research use only