

APIGENIN-6-GLUCOSIDE-8-ARABINOSIDE

Catalog No: #S3898


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	APIGENIN-6-GLUCOSIDE-8-ARABINOSIDE
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	564.496
Formulation	C26H28O14
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C1[C@@H]([C@@H]([C@H]([C@@H](O1)C2=C(C(=C(C3=C2OC(=CC3=O)C4=CC=C(C=C4)O)O)[C@H]5[C@@H]([C@H]([C@@H]([C@H](O5)CO)O)O)O)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in Pyridine, Methanol, EthanolAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Sun D,et al. Simultaneous determination of four flavonoids and one phenolic acid in rat plasma by LC-MS/MS and its application to a pharmacokinetic study after oral administration of the Herba Desmodii Styracifolii extract. J Chromatogr B Analyt Technol Biomed Life Sci. 2013 Aug 1;932:66-73.

Note: This product is for in vitro research use only