

Epimedin A

Catalog No: #S3906


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Epimedin A
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	838.8
Formulation	C39H50O20
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:CC1C(C(C(C(O1)OC2=C(OC3=C(C2=O)C(=CC(=C3CC=C(C)C)OC4C(C(C(C(O4)CO)O)O)O)C5=CC=C(C=C5)OC)O)OC6C(C(C(C(O6)CO)O)O)O)OPPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Pyridine, Methanol, Ethanol, etc.Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Li HB1, Chen F.Separation and purification of epimedin A, B, C, and icariin from the medicinal herb Epimedium brevicornum maxim by dual-mode HSCCC.J Chromatogr Sci. 2009 May-Jun;47(5):337-40.

Note: This product is for in vitro research use only