

Anemoside A3

Catalog No: #S3911


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anemoside A3
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	750.96
Formulation	C41H66O12
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OthersSMILES:C[C@H]1[C@@H]([C@H]([C@H]([C@@H](O1)O[C@@H]2[C@H]([C@H](CO[C@H]2O[C@H]3CC[C@@]4([C@H]5CC[C@@H]6[C@H]7[C@@H](CC[C@@]7(CC[C@@]6([C@@]5(CC[C@H]4[C@]3(C)CO)C)C(=O)O)C(=C)C)O)O)O)OPathways:OthersReceptor:OthersBo
ling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1.Ip FC, Fu WY, Cheng EY, Tong EP, et al.Anemoside A3 Enhances Cognition through the Regulation of Synaptic Function and Neuroprotection.Neuropsychopharmacology. 2015 Jul;40(8):1877-87.

Note: This product is for in vitro research use only