

## Complanatoside A

Catalog No: #S3922


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Complanatoside A
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	642.5
Formulation	C27H30O18
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:C1=C(C=C(C(=C1O)O[C@H]2[C@@H]([C@H]([C@@H]([C@H](O2)CO)O)O)O)C3=C(C(=O)C4=C(C=C(C=C4O3)O)O)O[C@H]5[C@@H]([C@H]([C@@H]([C@H](O5)CO)O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Pyridine, Methanol, Ethanol, etc.  
Appearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Zhang JJ1, Yan XL, Zhang YJ, et al.[Determination of complanatoside A in semen Astragali complanati by HPLC].Zhongguo Zhong Yao Za Zhi. 2005 Apr;30(8):600-2.

Note: This product is for in vitro research use only