

Astragaloside I

Catalog No: #S3925


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Astragaloside I
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	840.99
Formulation	C43H68O16
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others
 SMILES: CC(=O)O[C@H]1[C@@H](CO[C@H]([C@@H]1OC(=O)C)O[C@H]2CC[C@]34C[C@@]35CC[C@@]6([C@H]([C@H](C[C@]6([C@@H]5C[C@@H](C4C2)O)[C@H]7[C@@H]([C@H]([C@@H]([C@H](O7)CO)O)O)C)O)[C@]8(CC[C@H](O8)C(C)(C)O)C)O
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 10 mg/mL
 Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

- JIN Yi-qun, et al. Determination of astragaloside-I content in Qidan Jiannao Granules with HPLC-ELSD. Acta Chinese Medicine and Pharmacology, 2006-02

Note: This product is for in vitro research use only