

## Methyl L-histidinate dihydrochloride

Catalog No: #S3937


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Methyl L-histidinate dihydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	HDAC inhibitor
Calculated MW	242.1
Formulation	C7H13Cl2N3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: COC(=O)[C@H](CC1=CN=CN1)N.Cl. CIP pathways: Chromatin/Epigenetic; DNA Damage/DNA Repair;  
NF- $\kappa$ B Receptor: HDAC Boling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only