

L(+)-Homoarginine hydrochloride

Catalog No: #S3938


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	L(+)-Homoarginine hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	224.69
Formulation	C7H17CIN4O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:C(CCN=C(N)N)C[C@@H](C(=O)O)N.CIPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:water: 50 mg/mLAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only